



PROJECT REVIEW - 2

Artificial Intelligence-Enabled Surveillance of Knee Osteoarthritis Through Gait Analysis

Project Category: Research

Guide:

Dr. S. Karthick

Student Name & Registration Number:

Mr. Karthik B – RA2112701010014

Mr. Udhaya M – RA2112701010021

OBJECTIVE



Develop an AI-Based Gait Analysis Model

Implement computer vision and deep learning algorithms to extract movement features from video data.



Detect Post-Surgical Abnormalities and Deviations

Identify gait deviations, improper joint mobility, weight distribution issues, and other complications.



Automate the Post-Operative Assessment Process

Generate reports with key recovery metrics



DATASET OVERVIEW

- The KOA-NM dataset is designed for Knee Osteoarthritis (KOA) detection using gait analysis.
- It contains video recordings of individuals categorized based on KOA severity levels.
- A set of 6 red-colored passive reflective markers has been attached to the subject's body joints.
- The dataset is collected using a single NIKON DSLR 5300 camera placed 8m away from the walking mat.
- Videos are in .MOV format
- KOA-NM Dataset consists of gait videos classified into four categories:
 - NM (Normal Movement)
 - KOA_EL (Easy KOA)
 - KOA_MD (Moderate KOA)
 - KOA_SV (Severe KOA)



PRE PROCESSING

1. Video Format Conversion

- Converted .MOV files to .MP4 for compatibility.

2. Frame Extraction

- Extracted frames from videos at a uniform **FPS (Frames Per Second)**.

3. Noise Reduction & Enhancement

- Applied Gaussian Blur and contrast enhancement for better feature extraction.

4. Pose Estimation

- Used **Mediapipe** to extract joint movement key-points.

5. Key-Points Detection

- Identify the marker placed in joints of subject and store.

6. Data Augmentation

- Introduced rotation, flipping, and brightness variations to improve model generalization.



PRE PROCESSING

1. Challenges:

- Unreadable .MOV files
- Inconsistent FPS
- Blurred frames
- Missing joint key-points

2. Solutions:

1. Converted to .MP4
2. Standardized to fixed FPS
3. Applied sharpening & noise reduction
4. Improved pose estimation accuracy



FEATURE EXTRACTION

Input Data:

MOV format videos of patients walking **from left to right (side view)** across an **8-meter distance** in front of a fixed camera.

Pose Estimation with MediaPipe:

Extracts key body landmarks per frame:(x, y, z, visibility) for each joint.

Emphasis on lower-body joints (hip, knee, ankle) to assess post-knee replacement recovery.



FEATURE EXTRACTION

Preprocessing Steps:

Convert video to frames (e.g., 30 fps). Normalize and resize frames for uniformity.
Side view improves joint movement visibility during gait cycle.

Data Structuring:

Landmark data arranged as time series sequences:

Format: (frames, joints, features)

Example: 300 frames \times 33 joints \times 4 features

Padding/truncation applied for fixed sequence length.



Model Architecture

Architecture

1. Feature Extraction (CNNs) – Extracts spatial features from each frame.
2. LSTM Layer – Captures temporal dependencies in gait patterns.
3. Fully Connected Layer – Classifies severity levels of KOA.

LSTM for KOA Detection

- **Handles sequential gait data** effectively.
- **Remembers long-term dependencies** in walking patterns.
- **Outperforms traditional CNNs** in analysing movement over time.
- **Improves KOA severity classification** by learning time-series variations in gait

Architecture Diagram

